COLORED CLOSURES ONLY

Six styles of caps to choose from and two sealing types: O-ring seal and lip seal.

The cap is molded with a deep internal lip that fits snugly against the interior wall of the tube thus preventing the contents from coming in contact with the seal or threads, thereby reducing the chances of sample contamination. The cap's high profile facilitates manipulation especially in aseptic procedures and can remain attached (T332 & T336 Series) to the tube in order to eliminate mix-ups and contamination. All tubes and closures are manufactured in a clean environment.

Closures in natural color can be color.coded by the use of T345 Series colored Capinsert™ inserted on the top of the closure. This is accomplished without removing the cap. Colored caps are also available in all models as listed below.



Cat. #	Color	Qty/Pk
T340 NOS	Natural	1000
T340 BOS	Blue	1000
T340 GOS	Green	1000
T340 LOS	Lilac	1000
T340 ROS	Red	1000
T340 YOS	Yellow	1000
T340 WOS	White	1000





Cat. # Qty/Pk Color **T340 NLS** Natural 1000 Blue 1000 Green 1000

T340 BLS T340 GLS T340 LLS Lilac 1000 T340 RLS Red 1000 T340 YLS Yellow 1000 **T340 WLS** White 1000



Cat. #	Color	Qty/Pk
T340 NOSL	Natural	1000
T340 BOSL	Blue	1000
T340 GOSL	Green	1000
T340 LOSL	Lilac	1000
T340 ROSL	Red	1000
T340 YOSL	Yellow	1000
T340 WOSL	White	1000





Cat. #	Color	Qty/Pk
T340 NLSL	Natural	1000
T340 BLSL	Blue	1000
T340 GLSL	Green	1000
T340 LLSL	Lilac	1000
T340 RLSL	Red	1000
T340 YLSL	Yellow	1000
T340 WLSL	White	1000



Complete with o-ring seal and without loop. Made of polypropylene.

Cat. #	Color	Qty/Pk
T340 NOSFT	Natural	1000
T340 BOSFT	Blue	1000
T340 GOSFT	Green	1000
T340 LOSFT	Lilac	1000
T340 ROSFT	Red	1000
T340 YOSFT	Yellow	1000
T340 WOSFT	White	1000

The following closures have a flat top to accommodate automatic capping machines in packaging industry.





Lip seal and without loop.

Color	Qty/Pk
Natural	1000
Blue	1000
Green	1000
Lilac	1000
Red	1000
Yellow	1000
White	1000
	Natural Blue Green Lilac Red Yellow

T345

COLOR CODING CAPINSERT™ Made of polypropylene

Capinsert™ is used to color code a Micrewtube® according to your specific needs. They are inserted on the top of the closure and have a write-on frosted area for sample identification.



Cat. #	Color	Cat. #	Color	Qty/Pk	Qty/Cs
T345B	Blue	T345P	Pink	100	500
T345GY	Gray	T345R	Red	100	500
T345G	Green	T345V	Violet	100	500
T345L	Lilac	T345W	White	100	500
T345O	Orange	T345Y	Yellow	100	500
		T345AS	Assorted*	100	500

*Blue, lilac, red, yellow and white